


- MOUNT A5145, HOOK, AT 1372 [4'-6"] AFF. PROVIDE 2-A5145'S FOR NAVY PROJECTS ONLY. MOUNT 1-A5145 AT 1372 [4'-6"] AFF AND 1-A5145 AT 1524 [5'-0"] AFF.
- 2. COORDINATE LOCATION OF M1620, PATIENT CHART HOLDER, AND ROOM SIGNAGE.
 3. M4200, OTOSCOPE/OPHTHALMOSCOPE, IS OPTIONAL.
- 4. THE ACCEPTABLE RANGE FROM M5600, CHAIR, TO M5562, VISUAL ACCUITY TESTER, IS
- 3048 [10-0"] To 6096 [20"-0"].

 5. P3200, UITR ELECTRONIC LAVATORY, MAY BE USED INSTEAD OF P3100, LAVATORY. FOR THE P3200 PROVIDE EITHER BATTERY POWER OR HARD WIRED CONNECTION.
- 6. IF EQUIPMENT POLE IS MOUNTED ON THE M5016, REFRACTION DESK, INSTEAD OF BESIDE M5600, CHAIR, M5016 MAY NEED TO BE RELOCATED.

EYE LANE-FOLDED ELECTRONIC (EYEL4)

0 1000 2000 OPHTHALMOLOGY/OPTOMETRY FLOOR/EQUIPMENT PLAN 13.01 NSM (140 NSF) 1:50 SCALE (MILLIMETERS) DATE: 07/12/2001




KEY NOTES


1. SWITCH FOR ROOM ON/OFF LIGHTING CONTROLS. DIMMING CONTROLS ARE PROVIDED IN M5016, REFRACTION DESK.

EYE LANE-FOLDED ELECTRONIC (EYEL4)

1000 OPHTHALMOLOGY/OPTOMETRY 0 2000 CEILING PLAN 13.01 NSM (140 NSF) DATE: 07/12/2001 1:50 SCALE (MILLIMETERS)



NORTH



EAST

KEY NOTES

- KEY NOIES

 1. MOUNT A5145, HOOK, AT 1372 [4'-6'] AFF. PROVIDE 2-A5145'S FOR NAVY PROJECTS ONLY.

 MOUNT 1-A5145 AT 1372 [4'-6'] AFF AND 1-A5145 AT 1524 [5'-0"] AFF.

 2. P3200, UITR ELECTRONIC LAVATORY, MAY BE USED INSTEAD OF P3100, LAVATORY. FOR THE P3200 PROVIDE EITHER BATTERY POWER OR HARD WIRED CONNECTION.

 3. M4200, OTOSCOPE/OPHTHALMOSCOPE, IS OPTIONAL.

General Notes:

- 1. For typical equipment mounting heights see "TYPICAL EQUIPMENT MOUNTING HEIGHTS" drawing.
- 2. MIL—STD—1691 data can only be accessed from the Equipment Plan.

EYE LANE-FOLDED ELECTRONIC (EYEL4)

0 1000 2000 OPHTHALMOLOGY/OPTOMETRY **ELEVATIONS** 13.01 NSM (140 NSF) DATE: 07/12/2001 1:50 SCALE (MILLIMETERS)

As of 06/10/2005 Page 1

Room Code: EYEL4 Room Oty:	1 Room Area: 140.00 Room Description:	EYE LANE-FOLDED ELECTRONIC
NOOM COde. BIBLY NOOM OLV.	I NOOM ALEA. 140.00 NOOM DESCLIBETON.	DID DAND FOUND DUDCINONIC

JSN	NOMENCLATURE	Qty	UNIT ISSUE	LOG CAT	Utl 1	Utl 2	Utl 3	Utl 4	Utl 5	Utl 6			
Equipment per Room:													
A1066	Mirror, Float Glass, With SS Frame, 36X18	1	EA	А									
A5075	Dispenser, Soap, Disposable	1	EA	С									
A5080	Dispenser, Paper Towel, SS, Surface Mounted	1	EA	А									
A5106	Waste Disposal Unit, Sharps w/Glove Dispenser	1	EA	С									
A5145	Hook, Garment, Double, SS, Surface Mounted	1	EA	A					•				
F0205	Chair, Side With Arms	1	EA	С									
F2000	Basket, Wastepaper, Round, Metal, 18 H x 16 Dia.	1	EA	С									
F3200	Clock, Battery, 12" Diameter	1	EA	С					•				
M1620	Holder, Chart, Patient, Wall or Door Mounted	1	EA	A					•				
M4100	Sphygmomanometer, Aneroid, Wall Mounted	1	EA	С					•				
M4200	Otoscope/Ophthalmoscope, Wall Mounted	1	EA	С		A			•				
M5016	Desk, Refraction w/console, w/o Sink	1	EA	A		F							
M5030	Stool, Surgeon, Revolving	1	EA	С					•				
M5530	Lamp, Slit, w/Applanation Tonometer	1	EA	С		A							
M5535	Ophthalmoscope, Binocular, Indirect	1	EA	С		A							
M5545	Ophthalmometer/Keratometer	1	EA	С		A			•				
M5562	Tester, Visual Acuity, Video/Microprocessor Base	1	EA	С		A	•	•	•				
M5600	Unit, Exam, Eye, w/Motor Chair, Phoropter	1	EA	С		A							
M5700	Frame Set, Lens, Trial	1	EA	С				•	•				
P3100	Lavatory, Vitreous China, Slab Type	1	EA	A	D				•				
=======================================													

GRAND TOTAL